

Analytical Solutions

May 9, 2007 (Version 2.0)

for BioTechnology

BN 1214



Identification of Unknown Contaminants at a Subsurface Interface

Significance of Technique

The detection and identification of unknown, low-level contaminants present at an interface or in a buried layer is a challenging task that is often important to improve or troubleshoot semiconductor processes. These measurements require an analytical technique with the ability to detect a wide range of possible contaminants with high sensitivity.

For interfaces or layers in the near-surface region (i.e., the top 300 nm of the sample), these capabilities are provided by Time-of-Flight Secondary Ion Mass Spectrometry (TOF-SIMS). In these types of measurements, TOF-SIMS provides a survey analysis of either positive or negative ions with elemental detection limits as low as 10^8 atoms/cm².

In the example shown below a poly-silicon (30 nm)/polysilicon (30 nm) on silicon sample was analyzed by depth profiling. The cumulative mass spectrum acquired from the entire film structure indicated the presence of Cu, B, Na, K, Al and Ca. By reanalyzing the data file, the distributions of these elements could be reconstructed.



United States Locations

Tempe, Arizona
+1 480 239 0602 info.az@eaglabs.com
+1 602 470 2655 fax

Sunnyvale, California
810 Kifer Road
+1 408 530 3500 info.ca@eaglabs.com
+1 408 530 3501 fax

1135 E Arques Avenue
+1 408 738 3033
+1 408 530 3035 fax

785 Lucerne Drive
+1 408 737 3892
+1 408 737 3916 fax

Peabody, Massachusetts
+1 978 278 9500 info.ma@eaglabs.com
+1 978 278 9501 fax

Chanhassen, Minnesota
+1 952 828 6411 info.mn@eaglabs.com
+1 952 828 6449 fax

East Windsor, New Jersey
+1 609 371 4800 info.nj@eaglabs.com
+1 609 371 5666 fax

Syracuse, New York
+1 315 431 9900 info.ny@eaglabs.com
+1 315 431 9800 fax

Raleigh, North Carolina
+1 919 829 7041 info.nc@eaglabs.com
+1 919 829 5518 fax

Round Rock, Texas
+1 512 671 9500 info.tx@eaglabs.com
+1 512 671 9501 fax

International Locations

Shanghai, China
+ 86 21 6879 6088 info.cn@eaglabs.com
+ 86 21 6879 9086 fax

Tournefeuille, France
+ 33 5 61 73 15 29 info.fr@eaglabs.com
+ 33 5 61 73 15 67 fax

Frankfurt, Germany
+ 49 (0) 693053213 info.de@eaglabs.com
+ 49 (0) 69307941 fax

Tokyo, Japan
+ 81 3 5396 0531 info.jp@eaglabs.com
+ 81 3 5396 1930 fax

HsinChu, Taiwan
+ 886 3 5632303 info.tw@eaglabs.com
+ 886 3 5632306 fax

Uxbridge, United Kingdom
+ 44 (0) 1895 811194 info.uk@eaglabs.com
+ 44 (0) 1895 810350 fax